

His-tag Antibody-HRP labeled

Basic information:

Catalog No.:	UMA00163	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	ELISA: Use at an assay dependent dilution. WB: Use at an assay dependent dilution.
Reactivity:	N/A
Specificity:	Recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
Immunogen:	KLH-conjugated synthetic peptide HHHHHH.
Description:	Mouse monoclonal antibody to His-tag-HRP labeled
Uniprot:	
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data: