

B4GALNT1 Antibody


Basic information:

Catalog No.:	UPA00458	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human, Rat, Monkey
Specificity:	Recognizes endogenous levels of B4GALNT1 protein.
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human B4GALNT1. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to B4GALNT1
Uniprot:	Q00973
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of B4GALNT1 expression in Hela (A), C6 (B) whole cell lysates. (Predicted band size: 58 kD; Observed band size: 70 kD)