

CD274 Antibody


Basic information:

Catalog No.:	UPA02825	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB: 1:500 - 1:1000
Reactivity:	Human, Mouse, Rat
Specificity:	Recognizes endogenous levels of CD274 protein.
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CD274. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to CD274
Uniprot:	Q9NZQ7(Human), Q9EP73(Mouse)
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of CD274 expression in HEK293T (A), Raw264.7 (B), H9C2 (C) whole cell lysates.