

COX7B1 Antibody


Basic information:

Catalog No.:	UPA02870	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB: 1:500 - 1:1000
Reactivity:	Human, Mouse, Rat, Monkey
Specificity:	Recognizes endogenous levels of COX7B1 protein.
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human COX7B1. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to COX7B1
Uniprot:	P24311(Human), P56393(Mouse), P80431(Rat)
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of COX7B1 expression in HeLa (A), Raw264.7 (B), H9C2 (C) whole cell lysates.