

His-tag Antibody-HRP labeled


Basic information:

Catalog No.:	UPA05383	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:			

Useful Information:

Applications:	WB: 1:2000 - 1:5000
Reactivity:	N/A
Specificity:	Recognizes C-terminal, internal, and N-terminal His-tag fusion proteins.
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence of His-tag-HRP labeled. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to His-tag-HRP labeled
Uniprot:	
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of over-expressed His-tagged protein in 293T cell lysate.