

ECFP-tag Antibody


Basic information:

Catalog No.:	UPA05387	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB: 1:2000 - 1:5000
Reactivity:	N/A
Specificity:	Recognizes ECFP and CFP tag protein.
Immunogen:	Recombinant protein corresponding to ECFP-tag.
Description:	Rabbit polyclonal antibody to ECFP-tag
Uniprot:	
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of over-expressed ECFP-tagged protein in 293T cell lysate.