

Bovine Serum Albumin Antibody-HRP labeled

Basic information:

Catalog No.:	UPA05653	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	Immunogen affinity purified		

Useful Information:

Applications:	ELISA: Use at an assay dependent dilution.
Reactivity:	Bovine
Specificity:	Reacts only with Bovine Serum Albumin.
Immunogen:	Bovine Serum Albumin
Description:	Rabbit polyclonal antibody to Bovine Serum Albumin-HRP labeled
Uniprot:	
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data: