

CHD3 Rabbit mAb


Basic information:

Catalog No.:	UPA08555	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human
Specificity:	Recognizes endogenous levels of CHD3 protein
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within human CHD3. The exact sequence is proprietary.
Description:	Recombinant rabbit monoclonal antibody to CHD3
Uniprot:	Q12873
BiowMW:	Refer to Figures
Buffer:	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of CHD3 expression in HEK293T (A) whole cell lysates. (Predicted band size: 226 kD; Observed band size: 260 kD)