

ACTR1B Rabbit mAb


Basic information:

Catalog No.:	UPA08578	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human, Mouse, Rat
Specificity:	Recognizes endogenous levels of ACTR1B protein
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within human ACTR1B. The exact sequence is proprietary.
Description:	Recombinant rabbit monoclonal antibody to ACTR1B
Uniprot:	P42025
BiowMW:	Refer to Figures
Buffer:	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of ACTR1B expression in HEK293T (A), A549 (B), mouse kidney (C), rat kidney (D) whole cell lysates. (Predicted band size: 42 kD; Observed band size: 42 kD)