

ZNHIT6 Antibody


Basic information:

Catalog No.:	UPA09554	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human, Mouse, Rat
Specificity:	Recognizes endogenous levels of ZNHIT6 protein.
Immunogen:	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ZNHIT6. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to ZNHIT6
Uniprot:	Q9NWK9
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of ZNHIT6 expression in K562 (A), H1792 (B), mouse brain (C), mouse liver (D), rat brain (E) whole cell lysates. (Predicted band size: 53 kD; Observed band size: 54 kD)