

AK3 Antibody


Basic information:

Catalog No.:	UPA10198	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:2000)
Reactivity:	Human, Mouse, Rat
Specificity:	Recognizes endogenous levels of AK3 protein
Immunogen:	Recombinant fusion protein of human AK3. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to AK3
Uniprot:	Q9UIJ7
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of AK3 expression in A549 (A), mouse liver (B), mouse heart (C), rat pancreas (D) whole cell lysates. (Predicted band size: 18; 21; 25 kD; Observed band size: 26 kD)