

MSP58 Antibody


Basic information:

Catalog No.:	UPA10982	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:2000)
Reactivity:	Human, Mouse
Specificity:	Recognizes endogenous levels of MSP58 protein
Immunogen:	Recombinant fusion protein of human MSP58. The exact sequence is proprietary.
Description:	Rabbit polyclonal antibody to MSP58
Uniprot:	Q96EZ8
BiowMW:	Refer to Figures
Buffer:	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of MSP58 expression in MCF7 (A) whole cell lysates. (Predicted band size: 31; 50; 51; 53 kD; Observed band size: 55 kD)