

JAK2 Antibody


Basic information:

Catalog No.:	UPA13421	Source:	Rabbit
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human, Mouse
Specificity:	Recognizes endogenous levels of JAK2 protein.
Immunogen:	A synthetic peptide of human JAK2
Description:	Rabbit monoclonal antibody to JAK2
Uniprot:	O60674
BiowMW:	Refer to Figures
Buffer:	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Thimerosal and 0.05% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of JAK2 expression in K562 (A), NIH3T3 (B) whole cell lysates. (Predicted band size: 131 kD; Observed band size: 131 kD)