

PI3K p85 beta Antibody


Basic information:

Catalog No.:	UPA13625	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human, Mouse, Rat
Specificity:	Recognizes endogenous levels of PI3K p85 beta protein.
Immunogen:	Purified recombinant human PI3 Kinase p85 beta protein fragments expressed in E.coli.
Description:	Mouse monoclonal antibody to PI3K p85 beta
Uniprot:	O00459
BiowMW:	Refer to Figures
Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% Thimerosal, pH 7.3.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of PI3K p85 beta expression in C6 (A), K562 (B), MCF7 (C), U87MG (D), HCT116 (E), NIH3T3 (F), HeLa (G) whole cell lysates. (Predicted band size: 82 kD; Observed band size: 85 kD)