

DR5 Antibody


Basic information:

Catalog No.:	UPA13734	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was purified by immunogen affinity chromatography.		

Useful Information:

Applications:	WB (1:500 - 1:1000)
Reactivity:	Human
Specificity:	Recognizes endogenous levels of DR5 protein.
Immunogen:	Purified recombinant human DR5 protein fragments expressed in E.coli.
Description:	Mouse monoclonal antibody to DR5
Uniprot:	O14763
BiowMW:	Refer to Figures
Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% Thimerosal, pH 7.3.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of DR5 expression in SW480 (A), K562 (B), HCT116 (C), HeLa (D) whole cell lysates. (Predicted band size: 48 kD; Observed band size: 40,48 kD)