

GNPDA1 Antibody


Basic information:

Catalog No.:	UPA14352	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:			

Useful Information:

Applications:	WB (1:500 - 1:4000)
Reactivity:	Mouse
Specificity:	Recognizes endogenous levels of GNPDA1 protein.
Immunogen:	Recombinant fusion protein of human GNPDA1. The exact sequence is proprietary.
Description:	Mouse monoclonal antibody to GNPDA1
Uniprot:	
BiowMW:	Refer to Figures
Buffer:	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% Thimerosal.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.

Data:



Western blot analysis of GNPDA1 expression in mouse kidney (A) whole cell lysates. (Predicted band size: 32 kD; Observed band size: 35 kD)