

Anti-Hepatitis B Virus Core Antigen Antibody-FITC labeled

Basic information:

Catalog No.:	UDA00009	Source:	Goat
Size:	50ul/100ul	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	
Purification:	Immunogen affinity purified		

Useful Information:

Applications:	ELISA: Use at an assay dependent dilution.
Reactivity:	Hepatitis B virus
Specificity:	Reacts only with Hepatitis B Virus Core Antigen.
Immunogen:	Recombinant Hepatitis B Virus Core Antigen
Description:	Goat polyclonal antibody to Hepatitis B Virus Core Antigen-FITC labeled
Uniprot:	N/A
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage:	Store at 4°C short term and -28°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.