

Anti-Hepatitis B Virus Surface Antigen Antibody-HRP labeled

Basic information:

| | | | |
|-----------------------|-----------------------------|-------------------|------------|
| Catalog No.: | UDA00019 | Source: | Mouse |
| Size: | 50ul/100ul | Clonality: | Monoclonal |
| Concentration: | 1mg/ml | Isotype: | |
| Purification: | Immunogen affinity purified | | |

Useful Information:

| | |
|----------------------|--|
| Applications: | ELISA: Use at an assay dependent dilution. |
| Reactivity: | Hepatitis B virus |
| Specificity: | Reacts only with Hepatitis B Virus Surface Antigen. |
| Immunogen: | Recombinant Hepatitis B Virus Surface Antigen |
| Description: | Mouse monoclonal antibody to Hepatitis B Virus Surface Antigen-HRP labeled |
| Uniprot: | N/A |
| BiowMW: | N/A |
| Buffer: | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Storage: | Store at 4°C short term and -38°C long term. Avoid freeze-thaw cycles. |
| Note: | For research use only, not for use in diagnostic procedure. |