

# Anti-Hepatitis B Virus Core Antigen Anti-body-Biotin labeled

## Basic information:

<b>Catalog No.:</b>	UDA00023	<b>Source:</b>	Mouse
<b>Size:</b>	50ul/100ul	<b>Clonality:</b>	Monoclonal
<b>Concentration:</b>	1mg/ml	<b>Isotype:</b>	
<b>Purification:</b>	Immunogen affinity purified		

## Useful Information:

<b>Applications:</b>	ELISA: Use at an assay dependent dilution.
<b>Reactivity:</b>	Hepatitis B virus
<b>Specificity:</b>	Reacts only with Hepatitis B Virus Core Antigen.
<b>Immunogen:</b>	Recombinant Hepatitis B Virus Core Antigen
<b>Description:</b>	Mouse monoclonal antibody to Hepatitis B Virus Core Antigen-Biotin labeled
<b>Uniprot:</b>	N/A
<b>BiowMW:</b>	N/A
<b>Buffer:</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
<b>Storage:</b>	Store at 4°C short term and -42°C long term. Avoid freeze-thaw cycles.
<b>Note:</b>	For research use only, not for use in diagnostic procedure.