

Anti-Hepatitis B Virus Surface Antigen Antibody-Biotin labeled

Basic information:

Catalog No.:	UDA00025	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	
Purification:	Immunogen affinity purified		

Useful Information:

Applications:	ELISA: Use at an assay dependent dilution.
Reactivity:	Hepatitis B virus
Specificity:	Reacts only with Hepatitis B Virus Surface Antigen.
Immunogen:	Recombinant Hepatitis B Virus Surface Antigen
Description:	Mouse monoclonal antibody to Hepatitis B Virus Surface Antigen-Biotin labeled
Uniprot:	N/A
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage:	Store at 4°C short term and -44°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.