

Anti-Clenbuterol Antibody

Basic information:

Catalog No.:	UDA00028	Source:	Mouse
Size:	50ul/100ul	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	
Purification:	Protein G purified		

Useful Information:

Applications:	ELISA: Use at an assay dependent dilution.
Reactivity:	N/A
Specificity:	Reacts only with Clenbuterol.
Immunogen:	Clenbuterol conjugated to KLH
Description:	Mouse monoclonal antibody to Clenbuterol
Uniprot:	N/A
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage:	Store at 4°C short term and -47°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.