

## Anti-Deoxynivalenol Antibody

### Basic information:

<b>Catalog No.:</b>	UDA00032	<b>Source:</b>	Mouse
<b>Size:</b>	50ul/100ul	<b>Clonality:</b>	Monoclonal
<b>Concentration:</b>	1mg/ml	<b>Isotype:</b>	
<b>Purification:</b>	Protein G purified		

### Useful Information:

<b>Applications:</b>	ELISA: Use at an assay dependent dilution.
<b>Reactivity:</b>	N/A
<b>Specificity:</b>	Reacts only with Deoxynivalenol.
<b>Immunogen:</b>	Deoxynivalenol conjugated to KLH
<b>Description:</b>	Mouse monoclonal antibody to Deoxynivalenol
<b>Uniprot:</b>	N/A
<b>BiowMW:</b>	N/A
<b>Buffer:</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
<b>Storage:</b>	Store at 4°C short term and -51°C long term. Avoid freeze-thaw cycles.
<b>Note:</b>	For research use only, not for use in diagnostic procedure.