

## Rabbit Anti-Guinea Pig IgG (H&L) HRP

### Basic information:

<b>Catalog No.:</b>	US2124	<b>Source:</b>	Rabbit
<b>Size:</b>	100 ul/500 ul/1 ml	<b>Clonality:</b>	Polyclonal
<b>Concentration:</b>	1mg/ml	<b>Isotype:</b>	IgG
<b>Purification:</b>	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.		

### Useful Information:

<b>Applications:</b>	ELISA:1:5000-1:20000 WB:1:5000-1:20000 IHC:1:500-1:1000 ICC:1:500-1:1000
<b>Reactivity:</b>	Guinea Pig
<b>Specificity:</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
<b>Immunogen:</b>	Guinea Pig IgG
<b>Description:</b>	Rabbit Anti-Guinea Pig IgG (H&L) conjugated to HRP.
<b>Uniprot:</b>	N/A
<b>BiowMW:</b>	N/A
<b>Buffer:</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
<b>Storage:</b>	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
<b>Note:</b>	For research use only, not for use in diagnostic procedure.