

Mouse Anti-Monkey IgG (H&L) Biotin

Basic information:

Catalog No.:	US2239	Source:	Mouse
Size:	100 ul/500 ul/1 ml	Clonality:	Monoclonal
Concentration:	1mg/ml	Isotype:	
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.		

Useful Information:

Applications:	ELISA:1:10000-1:20000 WB:1:10000-1:20000 IHC:1:500-1:1000 ICC:1:500-1:1000
Reactivity:	Monkey
Specificity:	By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. Cross reactivity with human IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen:	Monkey IgG
Description:	Mouse Anti-Monkey IgG (H&L) conjugated to Biotin.
Uniprot:	N/A
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.