

Rabbit Anti-Bovine IgG (H&L) FITC

Basic information:

Catalog No.:	US2317	Source:	Rabbit
Size:	100 ul/500 ul/1 ml	Clonality:	Polyclonal
Concentration:	1mg/ml	Isotype:	IgG
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.		

Useful Information:

Applications:	IF:1:250-1:2000 FC:1:250-1:2000
Reactivity:	Bovine By immunoelectrophoresis and ELISA this antibody reacts specifically with
Specificity:	Bovine IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen:	Bovine IgG
Description:	Rabbit Anti-Bovine IgG (H&L) conjugated to FITC.
Uniprot:	N/A
BiowMW:	N/A
Buffer:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	For research use only, not for use in diagnostic procedure.